

Anti-Opri1 Polyclonal Antibody For Wb, Elisa - P41146

Item Number: CM0025835



Introduction

Rabbit polyclonal antibody against human OPR11 (nociceptin receptor, alias NOP, KOR-3) suitable for WB and ELISA applications. Useful in studying pain perception, nociception, locomotor activity, and G protein-coupled receptor signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	OPRL1 Rabbit pAb
Remarks/Alias	NOP, OOR, KOR3, NOPr, OPRL, ORL1, KOR-3, NOCIR, PNOCR, OPRL1
Species	Human
GeneID	4987
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 270-370 of human OPRL1 (NP_000904.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P41146
Protein Weight	41kDa
Shipping	Ice bag