

Anti-Opticin (Optc) Polyclonal Antibody For Wb And Elisa - Q9Ubm4

Item Number: CM0021353

Introduction

Rabbit polyclonal antibody against human OPTC (Opticin) suitable for WB and ELISA. Cross-reactivity with mouse and rat. Recombinant immunogen (a.a. 20-115). Ideal for angiogenesis and extracellular matrix studies.

[Learn More](#)

Parameter	Value
Product Name	OPTC Rabbit pAb
Remarks/Alias	OPT; OPTC
Species	Human
GenelD (Human)	26254
GenelD	26254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-115 of human OPTC (NP_055174.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UBM4
Protein Weight	37kDa
Shipping	Ice bag