


Anti-Optn Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q96Cv9

Item Number: CM0019170



Introduction

Rabbit polyclonal antibody against human OPTN (Optineurin), validated for Western blot, IHC-P, IP, and ELISA. Cross-reacts with mouse. Ideal for studying autophagy, innate immunity, and neurodegeneration.

[Learn More](#)

Parameter	Specification
Product Name	OPTN Rabbit pAb
Remarks/Alias	NRP; FIP2; HIP7; HYPL; ALS12; GLC1E; TFIIIA-INTP; OPTN
Species	Human
Gene ID (Human)	10133
Gene ID	10133
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-577 of human OPTN (NP_068815.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96CV9
Protein Weight	66kDa
Shipping	Ice bag