

Anti-Or10H3 Polyclonal Antibody For Wb, Ihc-P, Elisa - O60404

Item Number: CM0020324

Introduction

Rabbit polyclonal antibody targeting human olfactory receptor 10H3 (OR10H3), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for olfaction and GPCR research.

[Learn More](#)

Parameter	Value
Product Name	OR10H3 Rabbit pAb
Remarks/Alias	OR10H3
Species	Human
GeneID (human)	26532
GeneID	26532
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human OR10H3 (NP_039227.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60404
Protein Weight	36kDa
Shipping	Ice bag