

Anti-Or1A1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9P1Q5

Item Number: CM0020275

Introduction

OR1A1 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat olfactory receptor 1A1; ideal for GPCR and sensory transduction research.

[Learn More](#)

Parameter	Value
Product Name	OR1A1 Rabbit pAb
Remarks/Alias	OR17-7; OR1A1
Species	Human
GenelD (Human)	8383
GenelD	8383
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P1Q5
Protein Weight	35kDa
Shipping	Ice bag