

Anti-Or2T10 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Ngz9

Item Number: CM0014770

Introduction

Rabbit polyclonal antibody targeting human olfactory receptor OR2T10 (OR1-64), validated for Western blot and ELISA. Cross-reacts with mouse. Suitable for olfactory signaling and GPCR research studies.

[Learn More](#)

Parameter	Value
Product Name	OR2T10 Rabbit pAb
Remarks/Alias	OR1-64; OR2T10
Species	Human
GeneID (Human)	127069
GeneID	127069
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human OR2T10 (NP_001004693.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NGZ9
Protein Weight	35kDa
Shipping	Ice bag