

Anti-Orai1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96D31

Item Number: CM0025383



Introduction

ORAI1 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects ORAI1 (CRACM1) in human, mouse, rat. Ideal for studying store-operated calcium entry, immune signaling pathways, and channel function.

[Learn More](#)

Parameter	Value
Product Name	ORAI1 Rabbit pAb
Remarks/Alias	IMD9; TAM2; ORAT1; CRACM1; TMEM142A; ORAI1
Species	Human
Gene ID (Human)	84876
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ORAI1 (NP_116179.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96D31
Protein Weight	33kDa
Shipping	Ice bag