

Anti-Orc1 Polyclonal Antibody For Wb, Elisa - Q13415

Item Number: CM0016642



Introduction

High-quality rabbit polyclonal antibody against human ORC1, validated for WB and ELISA. Ideal for studying origin recognition complex function in DNA replication and pre-replication complex assembly.

[Learn More](#)

Parameter	Value
Product Name	ORC1 Rabbit pAb
Remarks/Alias	ORC1L; PARC1; HSORC1; ORC1
Species	Human
Gene ID	4998
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 696-856 of human ORC1 (NP_001177748.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13415
Protein Weight	97kDa
Shipping	Ice bag