


Anti-Orc6 Polyclonal Antibody For Wb, Elisa - Q9Y5N6

Item Number: CM0029560



Introduction

Rabbit polyclonal antibody targeting human ORC6, validated for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for DNA replication and origin recognition complex research.

[Learn More](#)

Parameter	Value
Product Name	ORC6 Rabbit pAb
Remarks/Alias	ORC6L; ORC6
Species	Human
Gene ID (Human)	23594
Gene ID	23594
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-252 of human ORC6 (NP_055136.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5N6
Protein Weight	28kDa
Shipping	Ice bag