


Anti-Orc6 Polyclonal Antibody (Rabbit) For Wb, Ihc-P, Elisa - Q9Y5N6

Item Number: CM0029267



Introduction

ORC6 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets origin recognition complex subunit 6 (28 kDa, Q9Y5N6), involved in DNA replication initiation. Ideal for nucleolar and DNA replication studies.

[Learn More](#)

Parameter	Value
Product Name	ORC6 Rabbit pAb
Remarks/Alias	ORC6L; ORC6
Species	Human
GeneID (human)	23594
GeneID	23594
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-252 of human ORC6 (NP_055136.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5N6
Protein Weight	28kDa
Shipping	Ice bag