

Anti-Orexin B/Hypocretin-2 Rabbit Polyclonal Antibody For Wb, Elisa - O43612

Item Number: CM0019079



Introduction

High-quality rabbit polyclonal antibody against human Orexin B/Hypocretin-2, validated for WB and ELISA. Ideal for studying feeding behavior, energy metabolism, and sleep-wake regulation. - O43612

[Learn More](#)

Parameter	Value
Product Name	Orexin B/Hypocretin-2 Rabbit pAb
Remarks/Alias	OX; PPOX; NRCLP1
Species	Human
GeneID (Human)	3060
GeneID	3060
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Orexin B/Hypocretin-2
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43612
Protein Weight	13kDa
Shipping	Ice bag