


Anti-Orm1/2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P02763 / P19652

Item Number: CM0002323



Introduction

Rabbit monoclonal antibody targeting human ORM1/2 (Alpha-1-acid glycoprotein 1/2), suitable for WB, IHC-P, and ELISA. Ideal for studying acute phase proteins, inflammation, and drug transport.

[Learn More](#)

Parameter	Value
Product Name	ORM1/2 Rabbit mAb
Remarks/Alias	AGP-B; AGP-B'; AGP2; ORM1/2
Species	Human
GeneID (Human)	5004/5005
GeneID	5004/5005
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P02763/P19652
Protein Weight	47kDa
Shipping	Ice bag