


Anti-Ormdl3 Rabbit Polyclonal Antibody For Wb And Elisa - Q8N138

Item Number: CM0017518



Introduction

This polyclonal antibody targets ORMDL3, an ER membrane protein essential for sphingolipid biosynthesis regulation. Validated for Western blot and ELISA, it reacts with human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	ORMDL3 Rabbit pAb
Remarks/Alias	ORMDL3
Species	Human
Gene ID (Human)	94103
Gene ID	94103
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ORMDL3 (NP_644809.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8N138
Protein Weight	17kDa
Shipping	Ice bag