


Anti-Orp1 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Bxw6

Item Number: CM0011974



Introduction

Rabbit monoclonal antibody targeting human ORP1 (OSBPL1A), suitable for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. ORP1 is a late endosomal sterol-binding protein involved in cholesterol transport and RAB7A regulation, ideal for lipid metabolism research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	ORP1 Rabbit mAb
Remarks/Alias	ORP1; ORP-1; OSBPL1B
Species	Human
GeneID (human)	114876
GeneID	114876
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 851-950 of human ORP1 (Q9BXW6).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXW6
Protein Weight	108kDa
Shipping	Ice bag