

Anti-Orp150/Grp170/Hsp12A Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y4L1

Item Number: CM0003180

Introduction

Rabbit monoclonal antibody against ORP150/GRP170/HSP12A (HYOU1) validated for WB, IHC-P, ELISA in human, mouse, rat. Suitable for endoplasmic reticulum stress and hypoxia research.

[Learn More](#)

Parameter	Specification
Product Name	ORP150/GRP170/HSP12A Rabbit mAb
Remarks/Alias	IMD59; Grp170; HSP12A; ORP150; GRP-170; ORP-150; ORP150/GRP170/HSP12A
Species	Human
GeneID (human)	10525
GeneID	10525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human ORP150/GRP170/HSP12A (Q9Y4L1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4L1
Protein Weight	111kDa
Shipping	Ice bag