

Anti-Osaldh2B1 Polyclonal Antibody For Wb, Elisa - Q9FrX7

Item Number: CM0020159

Introduction

OsALDH2B1 Rabbit Polyclonal Antibody (59 kDa) validated for WB and ELISA, specifically recognizing *Oryza sativa* ALDH2B1. Ideal for plant aldehyde dehydrogenase research and agricultural studies.

[Learn More](#)

Parameter	Value
Product Name	OsALDH2B1 Rabbit pAb
Remarks/Alias	Aldh2b; OsALDH2B1
Species	<i>Oryza sativa</i>
GeneID	4340725
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of <i>Oryza sativa</i> OsALDH2B1 (Q9FRX7)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q9FRX7
Protein Weight	59kDa
Shipping	Ice bag