

Anti-Osarf19 Polyclonal Antibody For Wb, Elisa - Q0D9R7

Item Number: CM0021804

Introduction

Rabbit polyclonal antibody targeting rice auxin response factor 19 (OsARF19, OsARF7a). Validated for Western blot (WB) and ELISA applications. Suitable for Oryza sativa research, involved in auxin signaling and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	OsARF19 Rabbit pAb
Remarks/Alias	ARF19; ARF7A; OsARF7a; OJ1215_E11.1; OsARF19
Species	Oryza sativa
GeneID	4341978
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 800-900 of Oryza sativa OsARF19 (Q0D9R7).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q0D9R7
Protein Weight	128kDa
Shipping	Ice bag