

## Anti-Osbip1 Polyclonal Antibody For Wb, Elisa - Q6Z7B0

Item Number: CM0022128

### Introduction

Rabbit polyclonal antibody targeting rice BiP1 (BIP1), an ER chaperone for protein folding and stress response. Validated for WB and ELISA in Oryza sativa. Ideal for seed development and ER stress research.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | OsBip1 Rabbit pAb                                                                                                      |
| Remarks/Alias    | BIP; OsBiP1; Osj_05115; OJ1442_E05.10; OsBip1                                                                          |
| Species          | Oryza sativa                                                                                                           |
| Gene ID          | 4328075                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of Oryza sativa OsBip1. (Q6Z7B0) |
| Source           | Rabbit                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                  |
| Application      | WB, ELISA                                                                                                              |
| Cross Reactivity | Oryza sativa                                                                                                           |
| SWISS            | Q6Z7B0                                                                                                                 |
| Protein Weight   | 73kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                |