


Anti-Osbpl2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H1P3

Item Number: CM0014886



Introduction

High-quality rabbit polyclonal antibody against human OSBPL2 (ORP2), validated for WB and ELISA. Ideal for studying lipid transport, cell adhesion, and deafness-related pathways.

[Learn More](#)

Parameter	Value
Product Name	OSBPL2 Rabbit pAb
Remarks/Alias	ORP2, ORP-2, DFNA67, DNFA67, OSBPL2
Species	Human
GeneID (Human)	9885
GeneID	9885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 291-480 of human OSBPL2 (NP_653081.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H1P3
Protein Weight	55kDa
Shipping	Ice bag