

Anti-Osbpl3 Polyclonal Antibody For Wb, Elisa - Q9H4L5

Item Number: CM0025023



Introduction

High-quality rabbit polyclonal antibody against human OSBPL3 (ORP3) for Western blot (WB) and ELISA applications. Recognizes human, mouse, rat OSBPL3. Recombinant immunogen from amino acids 370-520 of human OSBPL3.

[Learn More](#)

Parameter	Value
Product Name	OSBPL3 Rabbit pAb
Remarks/Alias	ORP3; ORP-3; OSBP3; OSBPL3
Species	Human
GeneID (Human)	26031
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 370-520 of human OSBPL3 (NP_056365.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4L5
Protein Weight	101kDa
Shipping	Ice bag