

Anti-Osbzip23 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Z312

Item Number: CM0022159



Introduction

Rabbit polyclonal antibody against rice OsbZIP23 transcription factor. Suitable for Western blot and ELISA, recognizing Oryza sativa target. Ideal for abiotic stress signaling and ABA pathway research.

[Learn More](#)

Parameter	Value
Product Name	Os bZIP23 Rabbit pAb
Remarks/Alias	ABF2; AREB1; BZIP23; OsABF2; OsAREB1; OsBZIP23; OJ1004_A11.20; OsbZIP23
Species	Oryza sativa
GeneID	4330838
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of Oryza sativa OsbZIP23. (Q6Z312).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q6Z312
Protein Weight	38kDa
Shipping	Ice bag