

Anti-Osein2 Polyclonal Antibody For Wb, Elisa - Q0D8I9

Item Number: CM0020799

Introduction

CamelBio supplies anti-OsEIN2 rabbit polyclonal antibody (Q0D8I9) reactive with *Oryza sativa*. Validated for WB and ELISA, targeting the ethylene signaling protein involved in rice development, senescence, and grain size.

[Learn More](#)

Parameter	Value
Product Name	OsEIN2 Rabbit pAb
Remarks/Alias	EIN2; MHZ7; OsEIN2; P0455F03.4
Species	<i>Oryza sativa</i>
Gene ID	4342431
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 576-710 of <i>Oryza sativa</i> OsEIN2 (Q0D8I9).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q0D8I9
Protein Weight	138kDa
Shipping	Ice bag