

Anti-Oself3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Snq6

Item Number: CM0023214

Introduction

Rabbit polyclonal antibody targeting OsELF3 (Q9SNQ6) from Oryza sativa. Validated for Western blotting and ELISA, suitable for flowering time regulation and photoperiodic pathway studies in rice.

[Learn More](#)

Parameter	Value
Product Name	OsELF3 Rabbit pAb
Remarks/Alias	Ef7; HD3B; Hd-q; Hd17; ELF3-1; OsELF3; OsELF3-1; OsJ_20086; OsELF3
Species	Oryza sativa
GeneID	4340087
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of Oryza sativa OsELF3. (Q9SNQ6).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q9SNQ6
Protein Weight	82kDa
Shipping	Ice bag