

Anti-Osfd1 Rabbit Polyclonal Antibody For Wb, Elisa - A0A0P0XR34

Item Number: CM0020333

Introduction

Rabbit polyclonal antibody against rice OsFD1 (P0478E02.1) for Western blot and ELISA. Targets a 21 kDa nuclear protein involved in abscisic acid signaling and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	OsFD1 Rabbit pAb
Remarks/Alias	P0478E02.1; OsFD1
Species	Oryza sativa
GeneID	9271661
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of Oryza sativa OsFD1. (A0A0P0XR34)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	A0A0P0XR34
Protein Weight	21kDa
Shipping	Ice bag