

Anti-Osfd2 Polyclonal Antibody For Wb And Elisa - Q5Z4G3

Item Number: CM0019417

Introduction

Rabbit polyclonal antibody against Oryza sativa OsFD2 (OsbZIP55). Validated for WB and ELISA. Suitable for studies in abscisic acid signaling, transcription regulation, and plant stress responses. Detects an 18 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	OsFD2 Rabbit pAb
Remarks/Alias	OsFD2; OsbZIP55
Species	Oryza sativa
Gene ID	9266154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of Oryza sativa OsFD2 (Q5Z4G3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q5Z4G3
Protein Weight	18kDa
Shipping	Ice bag