

Anti-Osgepl1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H4B0

Item Number: CM0028216

Introduction

High-quality rabbit polyclonal antibody targeting human OSGEPL1 (Q9H4B0) for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat, ideal for mitochondrial tRNA modification and genome maintenance studies.

[Learn More](#)

Parameter	Value
Product Name	OSGEPL1 Rabbit pAb
Remarks/Alias	Qri7; OSGEPL; OSGEPL1
Species	Human
Gene ID	64172
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human OSGEPL1 (NP_071748.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4B0
Protein Weight	45kDa
Shipping	Ice bag