

Anti-Osgin1 Polyclonal Antibody For Wb, Elisa - Q9Ujx0

Item Number: CM0013838

Introduction

Anti-OSGIN1 rabbit polyclonal antibody validated for Western blot and ELISA, targeting human oxidative stress-induced growth inhibitor 1 (OSGIN1), involved in cytokinesis and cell death regulation. Cross-reacts with human species. 52 kDa.

[Learn More](#)

Parameter	Value
Product Name	OSGIN1 Rabbit pAb
Remarks/Alias	BDGI; OKL38; OSGIN1
Species	Human
GeneID	29948
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human OSGIN1 (NP_892026.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UJX0
Protein Weight	52kDa
Shipping	Ice bag