

Anti-Osire1 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Xit1

Item Number: CM0020633

Introduction

OsIRE1 Rabbit Polyclonal Antibody validated for WB and ELISA in *Oryza sativa*. Targets serine/threonine-protein kinase/endoribonuclease IRE1 (99kDa), a key ER stress sensor. Suitable for studying unfolded protein response in rice.

[Learn More](#)

Parameter	Value
Product Name	OsIRE1 Rabbit pAb
Remarks/Alias	IRE1; OsIRE1; Osire1p; Osj_24190; OJ1103_E04.109
Species	<i>Oryza sativa</i>
Gene ID	4343198
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 436-650 of <i>Oryza sativa</i> OsIRE1 (Q7XIT1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q7XIT1
Protein Weight	99kDa
Shipping	Ice bag