

Anti-Osmads1 Rabbit Polyclonal Antibody For Wb, Elisa - Q10Pz9

Item Number: CM0022700

Introduction

Rabbit polyclonal antibody targeting *Oryza sativa* OsMADS1, validated for WB and ELISA. Suitable for investigating floral organ development, transcription factor function, and MADS-box protein research in rice.

[Learn More](#)

Parameter	Value
Product Name	OsMADS1 Rabbit pAb
Remarks/Alias	Lhs1; MADS1; OsMADS1; Osj_09913
Species	Oryza sativa
GeneID	4332059
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-257 of Oryza sativa OsMADS1 (Q10PZ9)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q10PZ9
Protein Weight	30kDa
Shipping	Ice bag