

Anti-Osmads15 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Q9I2

Item Number: CM0024286

Introduction

Rabbit polyclonal antibody against Oryza sativa OsMADS15, a MADS-box transcription factor. Validated for WB and ELISA, cross-reacts with rice. Ideal for plant developmental biology and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	OsMADS15 Rabbit pAb
Remarks/Alias	RAP1A; MADS15; FDRMADS3; OsMADS15; RMADS215; OJ1567_G09.109-1
Species	Oryza sativa
GeneID	4342214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 167-267 of Oryza sativa OsMADS15 (Q6Q9I2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q6Q9I2
Protein Weight	30kDa
Shipping	Ice bag