

Anti-Osmkk4 Polyclonal Antibody For Wb, Elisa - Q6K4Q0

Item Number: CM0018973

Introduction

Rabbit polyclonal antibody against rice OsMKK4, validated for Western blot (WB) and ELISA. Suitable for studying plant defense signaling and kinase function in *Oryza sativa*.

[Learn More](#)

Parameter	Value
Product Name	OsMKK4 Rabbit pAb
Remarks/Alias	OJ1316_E06.24; OsMKK4
Species	Oryza sativa
Gene ID	4330957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of Oryza sativa OsMKK4 (Q6K4Q0)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q6K4Q0
Protein Weight	40kDa
Shipping	Ice bag