

Anti-Osmpk6 Polyclonal Antibody For Wb And Elisa - Q336X9

Item Number: CM0023862

Introduction

High-quality rabbit polyclonal antibody targeting rice mitogen-activated protein kinase 6 (OsMPK6, Q336X9) for Western blot and ELISA. Ideal for studying MAPK pathways in *Oryza sativa*.

[Learn More](#)

Parameter	Value
Product Name	OsMPK6 Rabbit pAb
Remarks/Alias	MPK6; OsMPK6
Species	Oryza sativa
GeneID	4349225
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q336X9
Protein Weight	43kDa
Shipping	Ice bag