

## Anti-Osmr Polyclonal Antibody For Wb, Elisa - Q99650

Item Number: CM0030950

### Introduction

Rabbit polyclonal antibody against human OSMR, validated for WB and ELISA. Reacts with human and mouse. Ideal for studying oncostatin M and IL-31 signaling in immunology and dermatology research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | OSMR Rabbit pAb                                                                                                   |
| Remarks/Alias    | OSMRB; PLCA1; IL-31RB; OSMRbeta; IL-31R-beta; OSMR                                                                |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 9180                                                                                                              |
| GeneID           | 9180                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human OSMR (NP_003990.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | Q99650                                                                                                            |
| Protein Weight   | 111kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |