

Anti-Osnpr1 Polyclonal Antibody For Wb And Elisa - Q9Fdy4

Item Number: CM0022074

Introduction

Rabbit polyclonal antibody targeting Oryza sativa OsNPR1 (NPR1 homolog 1), validated for Western blotting and ELISA. Suitable for plant defense and disease resistance research.

[Learn More](#)

Parameter	Value
Product Name	OsNPR1 Rabbit pAb
Remarks/Alias	NH1; PR2; NPR1; NPR2; OsNH1; OsNPR1
Species	Oryza sativa
GeneID	4327315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 475-582 of Oryza sativa OsNPR1. (Q9FDY4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q9FDY4
Protein Weight	64kDa
Shipping	Ice bag