

## Anti-Osr1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Tax0

Item Number: CM0021510

### Introduction

High-quality rabbit polyclonal antibody targeting human OSR1 (odd-skipped-related 1), validated for WB, IHC-P, ELISA. Reacts with human, mouse, and rat. Suitable for embryonic development and transcriptional regulation research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | OSR1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | ODD; OSR1                                                                                                         |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 130497                                                                                                            |
| Gene ID          | 130497                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 28-115 of human OSR1 (NP_660303.1). |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q8TAX0                                                                                                            |
| Protein Weight   | 30kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |