

Anti-Osr2 Polyclonal Antibody For Wb, Elisa - Q8N2R0

Item Number: CM0024603

Introduction

Validated rabbit polyclonal antibody against OSR2 for WB and ELISA; recognizes human, mouse, and rat OSR2. Ideal for studying tooth germ development and zinc-finger transcription factors.

[Learn More](#)

Parameter	Value
Product Name	OSR2 Rabbit pAb
Remarks/Alias	OSR2
Species	Human
Gene ID (Human)	116039
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-312 of human OSR2 (NP_001135934.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N2R0
Protein Weight	36kDa
Shipping	Ice bag