

Anti-Osrack1A Polyclonal Antibody For Wb And Elisa - P49027

Item Number: CM0022202

Introduction

OsRACK1A rabbit polyclonal antibody validated for western blot (WB) and ELISA in rice (*Oryza sativa*). Detects the ~36kDa RACK1z protein, a positive regulator of innate immunity and plant defense. Ideal for rice blast research.

[Learn More](#)

Parameter	Value
Product Name	OsRACK1A Rabbit pAb
Remarks/Alias	RWD; GPB-LR; RACK1A; OsRACK1A
Species	Oryza sativa
GeneID	4324115
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of Oryza sativa OsRACK1A. (P49027)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	P49027
Protein Weight	36kDa
Shipping	Ice bag