

Anti-Ostc Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nrp0

Item Number: CM0019786



Introduction

OSTC Rabbit pAb validates in Western blot (WB) and ELISA for human, mouse, rat samples. Targets oligosaccharyltransferase complex subunit OSTC, key for N-glycosylation research and ER biology studies.

[Learn More](#)

Parameter	Value
Product Name	OSTC Rabbit pAb
Remarks/Alias	DC2; OSTC
Species	Human
Gene ID (Human)	58505
Gene ID	58505
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human OSTC (NP_067050.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRP0
Protein Weight	17kDa
Shipping	Ice bag