

Anti-Oste Rabbit Polyclonal Antibody For Wb, Elisa - Q10Sh0

Item Number: CM0022639

Introduction

Oste (TAD1/OsCCS52A) rabbit polyclonal antibody validated for Western blot and ELISA in *Oryza sativa*. WD-repeat cell cycle protein, 56 kDa. Ideal for rice research.

[Learn More](#)

Parameter	Value
Product Name	OsTE Rabbit pAb
Remarks/Alias	TAD1, TE, OsCCS52A, OsTE
Species	<i>Oryza sativa</i>
GeneID	4331448
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of <i>Oryza sativa</i> OsTE (Q10SH0)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q10SH0
Protein Weight	56kDa
Shipping	Ice bag