


Anti-Osteocalcin (Bglap) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P02818

Item Number: CM0014982



Introduction

Rabbit polyclonal antibody against human Osteocalcin (BGLAP), suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Targets bone Gla protein involved in bone mineralization and hormonal regulation. Ideal for IVD raw material sourcing.

[Learn More](#)

Parameter	Specification
Product Name	Osteocalcin (BGLAP) Rabbit pAb
Remarks/Alias	OC; BGP; OCN; Osteocalcin (BGLAP)
Species	Human
GeneID (Human)	632
GeneID	632
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Osteocalcin (BGLAP) (NP_954642.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02818
Protein Weight	11kDa
Shipping	Ice bag