


Anti-Osteocalcin/Ocn Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P02818

Item Number: CM0029350



Introduction

Rabbit polyclonal antibody targeting human Osteocalcin (BGLAP), suitable for WB, IF/ICC, and ELISA. Shows cross-reactivity with human, mouse, and rat. Supports bone biology and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	Osteocalcin/OCN Rabbit pAb
Remarks/Alias	OC; BGP; OCN; Osteocalcin (BGLAP)
Species	Human
Gene ID (Human)	632
Gene ID	632
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02818
Protein Weight	11kDa
Shipping	Ice bag