

Anti-Osteoglycin (Ogn) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P20774

Item Number: CM0030002

Introduction

Rabbit polyclonal antibody against human Osteoglycin (OGN) for Western blot, immunohistochemistry, immunofluorescence, ELISA. Cross-reacts with human, mouse, rat. Detects a ~34 kDa protein. Applicable for bone biology research.

[Learn More](#)

Parameter	Value
Product Name	Osteoglycin (OGN) Rabbit pAb
Remarks/Alias	OG; OIF; SLRR3A; Osteoglycin (OGN)
Species	Human
Gene ID (Human)	4969
Gene ID	4969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-298 of human Osteoglycin (OGN) (NP_054776.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20774
Protein Weight	34kDa
Shipping	Ice bag