


Anti-Osteopontin Monoclonal Antibody For Wb, Elisa - P10451

Item Number: CM0002224



Introduction

Rabbit monoclonal antibody targeting human Osteopontin (OPN/SPP1). Suitable for Western Blot and ELISA. Recognizes human Osteopontin; 35 kDa. Ideal for bone biology and immune response research.

[Learn More](#)

Parameter	Value
Product Name	Osteopontin Rabbit mAb
Remarks/Alias	OPN; BNSP; BSPI; ETA-1; Osteopontin
Species	Human
Gene ID (Human)	6696
Gene ID	6696
Immunogen	Synthetic peptide corresponding to amino acids 200-300 of human Osteopontin (P10451).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10451
Protein Weight	35kDa
Shipping	Ice bag