


Anti-Osteopontin Polyclonal Antibody For Wb, Elisa - P10451

Item Number: CM0013354



Introduction

Anti-Osteopontin polyclonal antibody validated for Western blot and ELISA. Detects human, mouse, rat osteopontin; recognizes amino acids 17-314. Ideal for studying cell-matrix interactions and cytokine functions.

[Learn More](#)

Specification	Detail
Product Name	Osteopontin Rabbit pAb
Remarks/Alias	OPN; BNSP; BSPI; ETA-1; Osteopontin
Species	Human
Gene ID	6696
Immunogen	Recombinant protein corresponding to amino acids 17-314 of human Osteopontin (NP_001035147.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10451
Protein Weight	35kDa
Shipping	Ice bag