


Anti-Ostm1 Polyclonal Antibody For Wb, Elisa - Q86Wc4

Item Number: CM0016778



Introduction

Rabbit polyclonal antibody against human OSTM1, also known as chloride channel 7 beta subunit, validated for WB and ELISA with cross-reactivity to mouse and rat. Essential for osteoclast and melanocyte maturation and function.

[Learn More](#)

Parameter	Value
Product Name	OSTM1 Rabbit pAb
Remarks/Alias	GL, GIPN, OPTB5, HSPC019, OSTM1
Species	Human
GeneID (human)	28962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-284 of human OSTM1 (NP_054747.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86WC4
Protein Weight	37kDa
Shipping	Ice bag