

Anti-Ostt1 Rabbit Polyclonal Antibody For Wb, Elisa - Q10Kf0

Item Number: CM0019473

Introduction

Rabbit polyclonal antibody against rice OsTT1 (proteasome subunit alpha type-2) validated for WB and ELISA. Recognizes Oryza sativa target with a predicted molecular weight of ~26kDa.

[Learn More](#)

Parameter	Value
Product Name	OsTT1 Rabbit pAb
Remarks/Alias	PAB1; OsPAB1; OsJ_11088; OsTT1
Species	Oryza sativa
Gene ID	4333000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 135-235 of Oryza sativa OsTT1 (Q10KF0)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q10KF0
Protein Weight	26kDa
Shipping	Ice bag