


Anti-Otc Rabbit Polyclonal Antibody For Wb And Elisa - P00480

Item Number: CM0026832



Introduction

Rabbit anti-OTC polyclonal antibody validated for Western blotting and ELISA. Detects human, mouse, and rat OTC. Essential for studying urea cycle, ammonia detoxification, and mitochondrial metabolism.

[Learn More](#)

Parameter	Value
Product Name	OTC Rabbit pAb
Remarks/Alias	OCTD; OTC1; OTCD; OTCas; OTC
Species	Human
GeneID (human)	5009
GeneID	5009
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 215-354 of human OTC (NP_000522.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00480
Protein Weight	40kDa
Shipping	Ice bag