

Anti-Otof Polyclonal Antibody For Wb, Elisa - Q9Hc10

Item Number: CM0021827

Introduction

OTOF rabbit polyclonal antibody validated for WB and ELISA in human samples. Recombinant immunogen (aa 1-300) targets otoferlin, a calcium sensor critical for synaptic vesicle exocytosis at cochlear hair cells, ideal for hearing and deafness research.

[Learn More](#)

Parameter	Value
Product Name	OTOF Rabbit pAb
Remarks/Alias	AUNB1; DFNB6; DFNB9; NSRD9; FER1L2; OTOF
Species	Human
GeneID (human)	9381
GeneID	9381
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human OTOF (NP_919224.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HC10
Protein Weight	227kDa
Shipping	Ice bag