


Anti-Otub1 Polyclonal Antibody For Wb, Elisa - Q96Fw1

Item Number: CM0016323



Introduction

High-quality rabbit polyclonal antibody targeting human OTUB1, a deubiquitinating enzyme involved in protein turnover, DNA repair, and immune regulation. Validated for Western blot and ELISA, cross-reacts with mouse. Ideal for basic research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	OTUB1 Rabbit pAb
Remarks/Alias	OTB1; OTU1; HSPC263; OTUB1
Species	Human
GeneID (Human)	55611
GeneID	55611
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human OTUB1 (NP_060140.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96FW1
Protein Weight	31kDa
Shipping	Ice bag