

Anti-Otub1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Fw1

Item Number: CM0016270



Introduction

Rabbit polyclonal antibody targeting human OTUB1, a Lys-48-specific deubiquitinase involved in protein turnover, DNA repair, and immune regulation. Validated for WB, IF/ICC, ELISA in human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	OTUB1 Rabbit pAb
Remarks/Alias	OTB1, OTU1, HSPC263, OTUB1
Species	Human
GeneID (Human)	55611
GeneID	55611
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FW1
Protein Weight	31kDa
Shipping	Ice bag